

Dictaphone PowerMic II Configuration in Dragon Medical 10

Description


The Dictaphone PowerMic II configuration screen is available in Dragon Medical 10 when a PowerMic II source profile has been created. These are the configuration options that are available under Tools, Options, PowerMic II.

Application Options:

- Global
- Dictation Box

Programmable Buttons:

- Transcribe
- Tab Backward
- Tab Forward
- Dictate
- Rewind
- Fast Forward
- Stop/Play
- Custom Left
- Enter/Select
- Custom Right



Button Programming Options

- **<None>** (*deactivate button*)
- **Applications**
 - Show DragonBar menu
 - Show dictation box
- **Dictation**
 - Microphone on/off
 - Show correction window
 - Microphone sleep/on
 - Pres-to-talk
 - Train new word
 - Transcribe recording
 - Transfer Text
- **Navigation**
 - Skip backward
 - Skip forward
 - Press OK button
 - Press Cancel button
 - Tab forward
 - Select previous variable field
 - Select next variable field
 - Clear delimiting characters
- **Recognition Modes**
 - Force normal recognition recognition
 - Force dictation recognition
 - Force command recognition
 - Force numbers recognition
 - Force spell recognition
 - Hidden Mode on/off
- **Playback**
 - Playback
 - Fast Playback
- **User Commands**
 - *Command used during the session*
 - *Command used during the session ...*
 - Specify


Application Options / Global Defaults




Dictation Box Defaults



Button Programming Options




Applications Options




Dictation Options




Navigation Options




Recognition Modes Options



Playback Options




User Commands* Options



***User commands** provides a list of the most recently used commands to choose from. To manually enter a command, select “**Specify ...**”

Specify User Command



For more information, please contact:

EXAQ, Inc.

PO Box 1325 Citrus Heights, CA 95611-1325
 (916) 966-8313 – Fax (916) 244-0582 - www.exaq.com - support@exaq.com